

Nuview Headband Visor System

Instructions for Use

Fitting the Standoff

1. Remove the blue Standoff from the packaging.
2. Set the loupes and light in accordance with the manufacturer's instructions to achieve good visualisation.
3. Clean the headband in the areas where the standoff contacts the headband to ensure a dry grease free surface



4. Offer the standoff up to the head band, aligning the centre marker with the centre of the loupes bar as shown as a guide to the definitive position, prior to adhering the Velcro.



5. Remove the protective tape from one of the white Velcro strips leaving the Velcro in situ.
6. Place the standoff in the definitive position and apply pressure to the standoff to ensure adhesion of the Velcro.
7. Remove the second Velcro protective tape and apply pressure to the standoff to ensure adhesion.
8. The standoff can now be removed and repositioned on the Velcro if minor realignment is required.

Re-ordering details

Pack of 2 Headband Standoffs
Pack of 10 Headband Visors

PPE-HB2001
PPE-HB2022

Nuview Ltd, Vine House, Selsley Rd, North Woodchester, Glos GL5 5NN

T +44(0)1453 872266 **F** +44(0)1453 872288

E info@nuview-ltd.com **W** www.nuview.co




Nuview Headband Visor System

Instructions for Use

Fitting the Visor

1. Remove the Visor from the packaging.
2. Press the visor onto the oblong standoff tabs as shown, making sure the visor sits flat against the standoff.



The visors are designed to be replaced after use and are supplied in packs of 10. Should the visor become scratched, discard the visor and replace with a new one. If cleaning and disinfecting the visor this should be undertaken in compliance with the practice/institution approved infection control procedures.

The Standoff can be removed for cleaning and disinfecting.

For storage with a visor attached invert the headband on a surface when not in use.

The Nuview Headband Visor system has been tested effective to EN 166:2003 and is CE marked accordingly.



Re-ordering details

Pack of 2 Headband Standoffs
Pack of 10 Headband Visors

PPE-HB2001
PPE-HB2022

Nuview Ltd, Vine House, Selsley Rd, North Woodchester, Glos GL5 5NN

T +44(0)1453 872266 **F** +44(0)1453 872288
E info@nuview-ltd.com **W** www.nuview.co

